

Intra Uterine Contraception


All traditional copper IUDs are the most effective Emergency Contraceptives, inserted up to 5 days after unprotected intercourse. All copper IUDs fitted to a woman over the age of 40 years may be left in situ until after the menopause. Recommended replacement intervals shown refer to the use of IUDs in younger women.



M Shape

Load 375

Copy of the familiar 'M-Shaped' device with 375mm² copper, for uterine length over 7cm.

	M/U	No	2	Black	5 yrs	32.5mm	19.5mm	3.6mm	>7cm	IUD0073001	GEZ171 x 3	£7.95

Frameless

GyneFix Viz

Unique frame-less device with 330mm² copper as tubes suspended on a surgical suture, the knot of which is embedded into the myometrium – special training required for fitting. Suitable for all uterine lengths. Stainless steel ring near knot ensures X-Ray visibility.

	Str.	N/A	1	Blue	5 yrs	30mm	2.2mm	4.75mm	Any size.	IUD0041200	-	£95.00

Plain Copper T

UT 380 Standard

Copy of the familiar 'copper T' design, loaded by pulling on the threads to draw the arms back into the introducer, 380mm² copper on vertical stem only. For uterine lengths 6.5 to 9cm.

	T	No	2	Black	5 yrs	32mm	32mm	3.8mm	6.5–9cm	IUD0072001	GEZ165 x 3	£8.95

References:


1. Professor J Guillebaud, Reviews in Gynaecological Practice, 1, (1), 2001, 24-31
2. Jo Dennis & Naomi Hampton, The Journal of Family Planning & Reproductive Healthcare, 28, (2), 2002, 61-68
3. Sunanda Gupta & Rosemary Kirkman, The Obstetrician & Gynaecologist, 1, (4), 2002, 37-43
4. NICE Clinical Guideline 30, October 2005

Intra Uterine Contraception



UT 380 Short

As UT 380 Standard (above), but now even smaller. Provides 380mm² copper suitable for minimum uterine length of 5cm.


	T	No	2	Black	5 yrs	27mm	32mm	3.8mm	>5cm	IUD0072002	GEZ166 x 3	£8.95

Banded Copper T

Copper T 380 devices with additional bands of copper on the horizontal arms give the highest efficacy and longest life in the patient. That's why they are regarded as 'The Gold Standard' by the leading authorities¹²³ and meet the NICE guidelines for LARC⁴.


TT 380 Slimline

'Banded copper T' design – the 'gold standard', featuring copper on horizontal arms as well as on the vertical stem, 380mm² copper, with easy-loading system and no capsule. For minimum uterine length of 7cm.

	T	Yes	2	Black	10 yrs	36.2mm	31.6mm	4.75mm	7cm (min)	IUD0071001	GEZ005	£9.95

Mini TT 380 Slimline

The only banded copper T device designed for use in the nulliparous. Features copper sleeves on horizontal arms as well as on the vertical stem, providing 380mm² copper on a short stem and with narrower arms. For maximum uterine length of 7cm.

	T	Yes	2	Black	5 yrs	29.5mm	23.2mm	4.75mm	7cm (max)	IUD0074001	GEZ160	£9.95

Intra Uterine Contraception


- Shape
- Bandage?
- Threads (No.)
- Threads (Col.)
- Life in patient
- Dimension (L)
- Dimension (W)
- Insertion width
- Cavity size
- Durbin codes
- NHS codes
- Price

T-Safe 380 A

The T-Safe CU 380A is a rebuild of the former, no longer available "Gyne-T".

The classic T-shaped IUD is the most common used IUD-model in the world, recommended even by the WHO. It convinces through secure fit in the uterus, lowest Pearl-Index and a minimal expulsion rate. Before insertion, the IUD has to be preloaded, which means that the horizontal arms have to be bent into the tube. For easy loading a loading capsule is enclosed.

Optimal, secure fit into the Uterus: minimal expulsion rate, Ideal alternative to the TT 380 Slimline.

	T	Yes	2	Blue	10yrs	35.9mm	31.9mm	4.5mm	6.6–9cm	IUD0021100 x 1	GBB002	£9.95
	T	Yes	2	Blue	10yrs	35.9mm	31.9mm	4.5mm	6.6–9cm	IUD0011100 (x 10)	GBB001	£97.00

T-Safe 380 QL

Identical characteristics to the T-Safe CU380. The T-Safe CU 380A QL convinces through secure fit in the uterus and in consequence low expulsion rate because of optimal memory test results.

Before insertion, the IUD has to be preloaded by using the plunger. Due to the development of the insertion body, loading of the IUD is possible in seconds.

	T	Yes	2	Blue	10yrs	35.9mm	31.9mm	4.75mm	7cm (Min)	IUD0024101 x 1	GEZ161	£10.45
	T	Yes	2	Blue	10yrs	35.9mm	31.9mm	4.75mm	7cm (Min)	IUD0024110 (x 10)	GEZ162	£100.25

Flexi T + 380


The Flexi-T IUD has a very simple design and a straightforward insertion method, however, one that is unique to the Flexi-T. The Flexi-T is ready for insertion as there is no plunger or elaborate applicator – the insertion tube fits around the stem of the IUD

	T	No	2	Blue	5yrs	32mm	28mm	3.5mm	6.5–9cm	IUD0033100 x 1	GEZ163	£10.06

Flexi T + 300

The Flexi-T+ 300 is the bigger version of the Flexi-T 300 (32 mm long and 28 mm wide, wrapped with 300 mm² of copper wire) and ideal for women who have given birth to one or more children, since research shows that the uterus is somewhat larger after pregnancy. The IUD is easy to insert, and insertion is almost completely painless.

Women who have a history of spotting throughout their menstrual cycle may want to opt for this type.

	T	No	2	Blue	5 yrs	32mm	28mm	3.5mm	6.6–9cm	IUD0032100 x 1	GEZ164	£8.47
										IUD0031100 (x 10)		£94.34

Intra Uterine Contraception

- Shape
- Bandage?
- Threads (No.)
- Threads (Col.)
- Life in patient
- Dimension (L)
- Dimension (W)
- Insertor width
- Cavity size
- Durbin codes
- NHS codes
- Price

Cu-Safe T 300 Copy of Flexi T

Same design as a Flexi T 300, Easy pain free insertion and removal. Secure double nylon thread system.

	T	No	2	Blue	5 yrs	22.9mm	22.9mm	3.5mm	6.6–9cm	IUD0015001 x 1	GEZ169	£9.71
										IUD0015001 (x 10)	GEZ170	£81.99

7MED 380 NsTA Standard (Copper-Silver)

An alternative to the highly popular UT380 Standard, loaded by pulling on the thread to draw the arms back into the introducer. 380mm of copper with a small thread of silver running through the middle.

	T	No	2	Blue	5 yrs	32mm	32mm	3.8mm	6.6–9cm	IUD0073002 x 1	GEZ183	£8.95


7MED 380 NsHa Short (Copper-Silver)

An alternative to the UT380 Short (page 7) but smaller.

	T	No	2	Blue	5 yrs	28mm	28mm	3.8mm	>5cm	IUD0073003 x 1	GEZ177	£8.95


Neo Safe

This model has the same design as the “Nova T 380”. Its copper surface amounts to 380 mm² guaranteeing safety up to 5 years. Also this model requires preloading before insertion. In contrast to the T-Safe CU 380A, the horizontal arms are to be pulled into the tube.

	T	No	2	Black	5 yrs	31.9mm	31.8mm	3.7mm	6.5–9cm	IUD0023100 x 1	GEZ167	£12.73
										IUD0014100 (x 10)	GEZ168	£117.66

Neo Safe Mini

The Neo-Safe T CU 380 Mini has identical features like the Neo-Safe T CU 380 and differs only in the dimensions from the standard model. Due to the less stem length as well as less width of the side arms, this model is ideally suited for women with a uterus length < 7 cm.

	T	No	2	Black	5 yrs	30mm	24mm	3.7mm	7 cm Max	IUD0023101 x 1	GB004	£12.73