

Manual Vacuum Aspiration (MVA) for uterine evacuation under local anaesthesia



49

Ipas EasyGrip Cannulae 49

Ⓢ

For use with the Aspirators on page 23. Individually wrapped, colour-coded by size.

| | | |
|-----------------------------------|-------------|-------|
| 3mm (for endometrial biopsy only) | DISP0022003 | £0.91 |
| 4mm | DISP0022004 | £1.39 |
| 5mm | DISP0022005 | £1.39 |
| 6mm | DISP0022006 | £1.39 |
| 7mm | DISP0022007 | £1.39 |
| 8mm | DISP0022008 | £1.39 |
| 9mm | DISP0022009 | £1.39 |
| 10mm | DISP0022010 | £1.39 |
| 12mm | DISP0022012 | £1.39 |



50

Ipas Karman Flexible Cannulae 50

Ⓢ

For use with the Ipas Aspirators. Individually wrapped.

| | | |
|------|-------------|-------|
| 4mm | DISP0230004 | £0.86 |
| 5mm | DISP0230005 | £0.86 |
| 6mm | DISP0230006 | £0.86 |
| 7mm | DISP0230007 | £0.86 |
| 8mm | DISP0230008 | £0.86 |
| 9mm | DISP0230009 | £0.86 |
| 10mm | DISP0230010 | £0.86 |
| 12mm | DISP0230012 | £0.86 |

See Adapter Kit below



51

Ipas Adapter Kit 51

Required to use Karman Cannulae sizes 4, 5, 6, 7, 8, 9 & 10mm with MVA Plus Aspirator and Double Valve Aspirator.

| | | |
|---------------------------|-------------|-------|
| Pack of 5 different sizes | DISP0240001 | £0.79 |
|---------------------------|-------------|-------|



52

Endometrial Biopsy Kit 52

Ⓢ

Contains one Single Valve Aspirator and 15 x 3mm Cannulae, all sterile.

| | |
|-------------|--------|
| DISP0291001 | £21.26 |
|-------------|--------|



Sexual Health Supplies

M.V.A.

for incomplete miscarriage & T.O.P. with local anaesthetic

- 98% effective
- Suitable from 5 to 12 weeks gestation
- May be carried out in OPD setting, freeing theatre for other cases
- May be used with Local Anaesthesia
- Easier identification of gestational sac to confirm evacuation



New RCOG Evidence-based Clinical Guideline No.7 (November 2011) states: Recommendation 7.2 "...electric or manual vacuum aspiration may be used as both are effective and acceptable to women and clinicians." Recommendation 7.14 "Services should be able to provide surgical abortion without resort to general anaesthesia"

For further information or to request a sample,
please call **020 8869 6590** or email **sexualhealth@durbin.co.uk**

www.durbin.co.uk

Manual Vacuum Aspiration (MVA) for uterine evacuation under local anaesthesia

MVA for surgical management of miscarriage, fetal demise and T.O.P.

98% effective
 Suitable from 5 to 12 weeks gestation
 May be carried out in OPD setting, freeing theatre for other cases
 May be used with Local Anaesthesia
 Easier identification of gestational sac to confirm evacuation
 Economical – theatre time estimated at £3,000+ per hour



53

Ipas Manual Aspirator

Patented locking valve system allows creation of vacuum prior to beginning the procedure, and finger-tip control of vacuum during the procedure.

| | | |
|------------------------|--------------------|---------------|
| Single Valve* | DISP0022060 | £6.56 |
| Double Valve 53 | DISP0021060 | £11.44 |

Ipas Sterile Double Valve Aspirator 53

DISP0025060 £12.51

Ipas Sterile Single Valve Aspirator*

DISP0025600 £7.83

Uterine Evacuation Kit

Contains one Double Valve Aspirator and four EasyGrip cannulae, sizes 6mm, 7mm, 8mm & 9mm, all sterile.

DISP0260102 £43.75

Single use sterile MVA Kit 54

Contains 1 x sterile Double Valve Aspirator, 1 x plastic Speculum broad, 1 x Luer Tenaculum, 1 x Dental Syringe, 1 x Rampley Spongeholder 24cm, 1 x Hegar Dilator 3/4, 1 x Hegar Dilator 5/6, 1 x Hegar Dilator 7/8, 1 x kidney dish, 1 x 60ml Gallipot, 1 x pair Nitrile gloves medium, 5 x large cotton wool balls, 1 x crepe paper 30 x 50cm, 1 x yellow waste bag and tie.

DISP0016019 £25.74

Although using MVA under LA offers maximum benefit, some units have also demonstrated significant savings of time and money using MVA under GA in theatre. To arrange a presentation on the use of the MVA equipment in your unit, please contact us.

*will accept cannulae sizes 4-6mm only.



54